

Complex Simplicity: The "Easy Button" Approach to Measurement Science

Tom Linz

Department of Chemistry
Wayne State University



Chemistry






Biomedical Opportunities

Biopharmaceuticals



Clinical Diagnostics



Liposomal Nucleic
Acid Drugs

Protein Drugs

Protein
Biomarkers

miRNA
Biomarkers




Analysis Performance Criteria

- Broadly applicable
- Inexpensive
- Rapid
- User-friendly
- High sensitivity
- High selectivity
- Low sample volume





Gel Electrophoresis

- Determines the chemical purity of biological samples
 - Identify which analytes are present
- Slab gels have limited analytical performance
 - LOD and resolution are insufficient






Electrophoresis Principles




Gel Electrophoresis

- Determines the chemical purity of biological samples
 - Identify which analytes are present
- Slab gels have limited analytical performance
 - LOD and resolution are insufficient
- Slab gel analysis sacrifices biological information
 - SDS denatures proteins prior to mass measurement
 - Cannot analyze particles




Thermal Gels

- Thermally responsive Pluronic polymers
 - Polymer micelles pack more densely at warm temperatures
 - Undergo a phase change
 - Viscosity is tunable with temperature



Thermal Gels



- Thermally responsive Pluronic polymers
 - Polymer micelles pack more densely at warm temperatures
 - Undergo a phase change
 - Viscosity is tunable with temperature



Adapted from *J. Pharm. Pharm. Sci.*
2006, 9, 339-358


Microfluidic Device Operation

- Soft lithography used to fabricate PDMS microfluidic devices
 - 100 x 20 μm channels
- Sample is directly cast into thermal gel
- Images acquired using an epifluorescence microscope
- Temperature controller integrated into microscope stage




Microfluidic Separations

Gated Electrophoresis




Isotachopheresis



Thermal Gel Electrophoresis (TGE)



- Use simple single-channel microfluidic devices
 - Cast sample in thermal gel
 - Load thermal gel throughout channel
- Inline preconcentration and separation
 - No need to direct multiple fluid flows
 - No sample injection
 - No precise voltage switching
 - No user intervention






miRNA Analysis

- MicroRNAs silence genes
 - Prevent protein translation from mRNA
 - Misregulation causes pathogenesis
 - Biomarkers for diseases
 - Being developed into therapeutics


- miRNAs composed of ~ 22 nucleotides
 - Same-length sequences do not separate in gel electrophoresis
 - Cannot distinguish different species



Multiplexed miRNA TGE




Multiplexed miRNA TGE




LOD \approx 100 pM

Second Generation TGE Devices





LOD \approx 100 pM



LOD \approx 10 pM


Native Protein Analysis

- Proteins must remain correctly folded to exhibit proper bioactivity
 - Either inactive or harmful if misfolded or aggregated
- Proteins must have proper multimeric units or complexes associated
 - Post-translation modifications must be present



Native Protein TGE

- Fluorescent ovalbumin used as a model protein
 - Multiple folding and phosphorylation variants




Analysis Performance Criteria

- Broadly applicable ✓
- Inexpensive ✓
- Rapid ✓
- User-friendly ✓
- High sensitivity ✓
- High selectivity ✓
- Low sample volume ✓




Vesicle Analysis

- Vesicles are lipid nanoparticles that contain biomolecules
- Liposomes are drug delivery vehicles that protect internal cargo
- Extracellular vesicles are secreted by cells for communication
- Need to interrogate intravesicular biomolecules





10 °C




Joule Heating Barriers

?



Cell Enrichment




HEK 293 cells


Stage = 10 °C

10x speed


Cell Enrichment




HEK 293 cells
Stage = 10 °C
10x speed



Enrich





Collect



■ 10 °C
■ 22 °C

Conclusions

- Thermal gel electrophoresis provides inline enrichment and separation
 - Analyze proteins, nucleic acids, and vesicles in a user-friendly analysis
- Digital quantitative PCR provides single-molecule sensitivity and a large dynamic range
 - Quantify vesicles AND intravesicular biomolecules in a single analysis
 - Measure extracellular biomolecules with good LODs



Acknowledgements

- Linz Lab members
 - Bailey McCarthy Riley
 - Shakila Peli Thanthri
 - Mario Cornejo
 - Hao Mai
 - Courtney Cunningham
 - Brice Vanness



- Funding
 - NIH R21GM137278
 - NSF 2046487
 - SACP Starter Grant
 - Wayne State University




National Institute of
General Medical Sciences




Microfabrication


1. Mold Fabrication

- 1 SU-8 photoresist application to Si wafer




SU-8 photoresist
Si wafer
- 2 Spin-coating and soft baking


- 3 Photomask alignment and UV exposure




UV light
Photomask
- 4 SU-8 development, baking, and rinsing





2. PDMS Device Fabrication

- 5 PDMS casting and thermal annealing



70 °C
- 6 PDMS chip peeling-off


- 7 Creating ports


- 8 Lay glass slide onto PDMS chip

